New Safety Count System

Disarmer of scalpel blades and container of needles and blades

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Instructions for use

Preparing the product for use:

Open the sterile packaging according to standard sterile-field procedure and position the product either on the sterile table

Using the product:

- ◆ To store the needles used during surgery, the Operator, based on the model of container in use, can decide whether simply allow the needles to drop onto the adhesive foam (when present) (after having removed the protective film) (fig. A and A1); or place them on the foam bar by perforating it along the numbered "needle counter" scale (when present) (fig. B); or position them along the magnetic, numbered "needle counter" scale (when present) (fig. C).
- ◆ To detach used scalpel blades without needing to touch them, use the convenient two-sided blade detacher which is standard on all code numbers, following these specific instructions (fig. D):
 1. Hold the scalpel with mounted blade in the same manner as during the blade attachment procedure, with the attachment mechanism upwards.
 2. In this position, insert the blade into the blade detachment opening, stopping once the blade is entirely inserted.
 3. With the thumb, apply downward pressure at the thinnest point of the scalpel handle's neck, vertically pushing the non-cutting distal end of the blade until it lifts up from its attached position.
 4. While applying pressure as described in point 3, pull the scalpel handle back horizontally. The scalpel blade will detach and remain inside the protective collection housing.
- ◆ Once the surgery has been completed and, if necessary, having checked the correspondence of the number of needles and blades contained in the NEW SAFETY COUNT SYSTEM with that of the open envelopes on the sterile table, close the container as follows: lift one of the two halves and bring it to the other, checking that no sharps protrude from the perimeter. Then exert a slight pressure until you hear a "click" (fig.E), which guarantees the closure of the container, which can in any case be reopened for a further check of the contents.

How to reopen the SAFETY COUNT:

Using the scalpel handle, insert the tip into the detachment hole (fig. F), and while thumb and forefinger of the other hand hold the upper side firm, apply a light lever action in an outward direction until it becomes possible to lift the upper side.

Warnings

- ⇒ Do not use if the packaging is damaged. The product must be used immediately after opening the package.
- ⇒ Product for surgical use. To be used only by adequately trained personnel.
- ⇒ Do not fill beyond the maximum capacity: all blades and needles must be positioned correctly so that they are secured. In particular, the needles must always be inserted flat into the container, ie parallel to the horizontal surface.
- ⇒ None of the sharps must protrude over the edges of the device so as not to hinder its closing.
- ⇒ Store in a dry place at room temperature.

All serious incidents occurring in relation to this medical device should be reported to the Manufacturer and the Competent Authority of the country where the user or patient is established, as required by European Regulation 745/2017.

Disposal

The device must be disposed of as hospital waste according to the laws in force.

Storage conditions

To maintain the correct conservation over time it is advisable to always use the original packaging together with a dry environment with temperatures between + 1 °C and + 40 °C, not in direct contact with heat sources, protected from dust and harmful materials.

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MED EUROPE S.r.l.