Sterile Doppler Probe

DESCRIPTION

Doppler probes specific for the vascular sectors, with various diameters, frequencies; 8, 10 or 20MHz, sterile, disposable, in double sterile packaging.

INTENDED USE

The device acts as a transducer during surgery operations to signal and detect the presence of blood flow inside the blood vessels.

To be used for a period of time up to 60 minutes maximum.

WARNINGS

Before using the device, carefully read the instructions for use. The cannula must be used exclusively in the hospital setting, inserted, handled and removed exclusively by a qualified medical doctor or under explicit instructions by the same.

The medical techniques described in these instructions are supplied purely as information. The surgeon must determine the adequacy of the procedure suggested, based on his/her own experience and the conditions of the patient. The surgical equipe is solely responsible for the management of the surgical techniques that are appropriate.

The probes must be used exclusively with the doppler models indicated below; any other type of connection could cause injury to the patient and/or the user.

This device is not intended, sold or to be utilised in any other way other than that indicated;

The product is single use disposable, therefore no adjustments or repairs are allowed on the instrument, so as to avoid injury to the patient and the user;

All serious incidents occurring in relation to this medical device should be reported to the Manufacturer and the Competent Authority of the country where the user or patient is established, as required by European Regulation 745/2017.

APPLICATION SPECIFICATIONS

Available in frequencies of 8, 10 or 20MHz, in diameters 0.8-1,2-2.0 mms.

Vascular flow detection depth based on frequency:
10mm./8MHz, 5mm./10MHz, 3mm./20MHz.

Compatible with HADECO Doppler models:
Smartdop SD-45 and Bidop7

CONTRAINDICATIONS

- Do not reuse or resterilise as this is a disposable product. Reuse of the device determines a high risk of bacterial contamination and cross infections;
- The device should not be used when in the opinion of the Doctor it may be detrimental to the patient's health.

PRECAUTIONS

- The contents of the packaging are sterile, if intact, undamaged and closed; if not intact the device must not be used;
- The product must be used immediately after opening the package;
- After opening the package, examine the product carefully: do not use if defects are found;
- When available, choose the frequency according to the typology of application.

INSTRUCTIONS FOR USE

- Select the probe model and its frequency based on the depth of the vessel under investigation;
- Use sterile aseptic procedure to remove the probe from the package;
- Connect the cable of the probe to the doppler in use;
- Verify the functioning of the connection to the doppler, checking the noise emission;
- Adhere and orient the tip of the probe on the tissues overlying the vessel under investigation, always making sure that the tip of the probe is moistened either with the patient's blood exudate or with sterile saline.
- For further indications, warnings and mode of function of the doppler, refer to the relative instruction Manual.

DISCLAIMER OF LIABILITY

We declare under our responsibility that the product complies with EU Regulation 745/2017 for medical devices.

MED EUROPE® is not held responsible for any type of loss, damage or injury of any kind, direct or indirect, or for damage resulting from the improper use of its products.

In particular, MED EUROPE® is not required to reimburse the costs incurred by the buyer or third parties following the improper use of the device. This also applies to any medical expenses, hospital treatment, medications or consequential damages.


DISPOSAL

The device must be disposed of as hospital waste according to the laws in force.

STORAGE CONDITIONS

To maintain the correct conservation over time it is advisable to always use the original packaging together with a dry environment with temperatures between +1 °C and +40 °C, not in direct contact with heat sources, protected from dust and harmful materials.

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MED EUROPE S.r.l.