


# Nasal Splint Vented



Following a rhinoseptoplasty surgery, patients can't breathe through their nose. Nasal Splint Vented allows a patient to breathe normally while providing effective support to the septum. It is made of 100% silicone to guarantee maximal biocompatibility and non-traumaticity.

- Made of "medical grade" silicone, Shore 60.
- Anatomical design, ideal for septoplasty and rhinoplasty procedures.
- Sterilised by gamma rays.
- Latex free.
- Packaged in double envelope.

REF MED 7370

2 pieces per envelope, 5 envelopes per box.

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
 www.simitalia.it www.medeurope.eu  
 E-Mail: info@medeurope.eu

Distributed by:

