


# Retractor for NEUROSURGERY

## Hook for Galea

Revision date: 26.01.2017




### Description

Hooked retractor for galea for use in neurosurgery.

### Characteristics

- Made of AISI 316L steel
- Diameter 1,6mm.
- "Medical grade" silicone elastic attachment
- Latex free device
- Packaged in double sterile envelope
- Gamma ray sterilization
- Box of 50 pcs.



### Codes

REF MEDHOOK001

Large size, 110mms.

REF MEDHOOK002


Large size, 110mms., coated with "medical grade" silicone

REF MEDHOOK003

Small size, 100mms.

REF MEDHOOK004

Small size, 100mms., coated with "medical grade" silicone



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
www.simitalia.it www.medeurope.eu  
E-Mail: info@medeurope.eu

Distributed by:



# Retractor for NEUROSURGERY

## Double Hook for Galea

Revision date: 22.10.2020



### Description



Double hook retractor for galea for use in neurosurgery.

### Characteristics

- Made of AISI 301 steel
- Thickness 0,7 mm
- Length 33 mm
- "Medical grade" silicone elastic attachment
- Latex free device
- Packaged in double sterile envelope
- Sterilised by gamma rays
- Box of 50pz.

### Code

REF MEDHOOK-005  
(double hook)



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
www.simitalia.it www.medeurope.eu  
E-Mail: info@medeurope.eu

Distributed by:

