

Disposable Doppler Probe

Flexible and malleable, sterile disposable model




CLASS IIa

0051

RULE 6

STERILE R



Flexible or malleable distal element
Length of 120 mms.


Ergonomic and slip-resistant handle
with medical grade silicone insert

Characteristics

- Available in two frequencies of 10 or 20 MHz
- Flexible or malleable distal element
- Ergonomic and slip-resistant handle with medical grade silicone insert
- Total length of 250 mms.
- Distal element's length of 120 mms.
- Tip diameter: 0,8 - 1,2 - 2 mms.
- Disposable
- Handle made of styrolux
- Sterilised by Gamma Rays
- Double sterile packaging
- 6 pcs. Box




Flexible and malleable probe, 20 MHz frequency




Optional: millimeter
scale at the tip

Malleable probe, 20 MHz frequency

Compatible with HADECO Doppler models:
SD-45, BIDOP 7



Flexible and malleable probe
Length 400 mm



Pencil probe, 8 MHz frequency

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



REF.	Frequency	Focal depth	Tip Diameter	Handle colour	Model
MED NDP-08-2.0*	8 MHZ	10 mm	2 mm	GREEN	Pencil probe
MED NDP-10-2.0*	10 MHZ	5 mm	2 mm	WHITE	
MED NDP-20-2.0*	20 MHZ	3 mm	2 mm	BLUE	
MED NDP-10-0.8-F*	10 MHZ	5 mm	0.8 mm	WHITE	Flexible and malleable probe
MED NDP-20-0.8-F	20 MHZ	3 mm	0.8 mm	BLUE	
MED NDP-08-1.2-F*	8 MHZ	10 mm	1.2 mm	GREEN	
MED NDP-10-1.2-F*	10 MHZ	5 mm	1.2 mm	WHITE	
MED NDP-20-1.2-F	20 MHZ	3 mm	1.2 mm	BLUE	
MED NDP-08-2.0-F*	8 MHZ	10 mm	2 mm	GREEN	
MED NDP-10-2.0-F*	10 MHZ	5 mm	2 mm	WHITE	
MED NDP-20-2.0-F	20 MHZ	3 mm	2 mm	BLUE	
MED NDP-08-2.0-Z*	8 MHZ	10 mm	2 mm	GREEN	
MED NDP-10-2.0-Z*	10 MHZ	5 mm	2 mm	WHITE	Flexible and malleable probe Distal length 180 mm
MED NDP-20-2.0-Z*	20 MHZ	3 mm	2 mm	BLUE	
MED NDP-10-1.2-M*	10 MHZ	5 mm	1.2 mm	WHITE	
MED NDP-20-1.2-M	20 MHZ	3 mm	1.2 mm	BLUE	Malleable probe
MED NDP-08-2.0-M*	8 MHZ	10 mm	2 mm	GREEN	
MED NDP-10-2.0-M*	10 MHZ	5 mm	2 mm	WHITE	
MED NDP-20-2.0-M	20 MHZ	3 mm	2 mm	BLUE	
MED NDP-08-2.0-L*	8 MHZ	10 mm	2 mm	GREEN	
MED NDP-10-2.0-L*	10 MHZ	5 mm	2 mm	WHITE	Malleable probe with millimeter scale at the tip
MED NDP-20-2.0-L	20 MHZ	3 mm	2 mm	BLUE	
MED NDP-08-2.0 LAP*	8 MHZ	10 mm	2 mm	GREEN	
MED NDP-10-2.0 LAP*	10 MHZ	5 mm	2 mm	WHITE	Flexible and malleable probe Distal length 400 mm
MED NDP-20-2.0 LAP*	20 MHZ	3 mm	2 mm	BLUE	

* Product not in stock; available by specific request