

# MED X-Ray Drain

## Drainage with radiopaque line




### Characteristics

Silicone drainage for use in the operating room. Produced in soft, biocompatible and hydrorepellent silicone.

Considering that the new European Directive 93/42 has imposed limitations on drainages created on site in the operating room which are made from PVC or Medical Grade Silicone in rolls, the need has been created for the medical industry to develop a series of special drainages, already sterile and ready for use, which are capable of satisfying the needs encountered in the operating room. These types of drainages are particularly used in abdominal and thoracic surgery, and have different requisites from standard drainages both with regard to length, softness and the extension of perforations.

### Specifications

Packaged singularly in double sterile (gamma ray) packaging. One box contains 20 pcs.



|  | Int. Ø<br>mm. | Ext. Ø<br>mm. | Number<br>of holes | Hole Ø<br>mm. | Shore | Ext. CH | REF<br>tube length 70<br>cm. | REF<br>tube length 50<br>cm. |
|--|---------------|---------------|--------------------|---------------|-------|---------|------------------------------|------------------------------|
|  | 3,35          | 5,5           | 8                  | 4             | 62    | 15      | A03350550C6T                 | A03350550C6R                 |
|  | 3,18          | 6,35          | 8                  | 4             | 62    | 19      | A03180635C6T                 | A03180635C6R                 |
|  | 5             | 7             | 8                  | 5             | 62    | 21      | A05000700C6T                 | A05000700C6R                 |
|  | 4,8           | 7,9           | 8                  | 5             | 62    | 24      | A04800790C6T                 | A04800790C6R                 |
|  | 6,4           | 9,5           | 8                  | 7             | 62    | 30      | A06400950C6T                 | A06400950C6R                 |
|  | 6,4           | 9,5           | 16                 | 3             | 62    | 30      | A06400950G6T*                | A06400950G6R                 |
|  | 6,4           | 11            | 8                  | 8             | 62    | 33      | A06401100C6T                 | A06401100C6R                 |
|  | 6,4           | 11            | 16                 | 3             | 62    | 33      | A06401100G6T*                | A06401100G6R                 |
|  | 9,5           | 12,7          | 8                  | 8             | 62    | 36      | A09501270C6T*                | A09501270C6R                 |
|  | 7,9           | 12,7          | 8                  | 8             | 62    | 36      | A07901270C6T                 | A07901270C6R                 |
|  | 7,9           | 12,7          | 16                 | 3             | 62    | 36      | A07901270G6T*                | A07901270G6R                 |

\* Product not in stock, available by specific request

Specific drainages may be produced upon request, depending on the quantity requested, using silicone tubes of various lengths, numerous CH sizes and with different diameters and number of holes.

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
www.simitalia.it www.medeurope.eu  
E-Mail: info@medeurope.eu

Distributed by:

